FR Informations et instructions concernant le produit / EN Product Information and Instructions / IT Informazioni sul prodotto ed istruzioni per l'uso / DE Produktinformationen und Gebrauchsanweisung / ES Instrucciones e información sobre el producto / NL Productinformatie en instructies / DA Produktinformation og -vejledning / FI Tuotetta koskevat tiedot ja käyttöohjeet / SV Produktinformation och instruktioner / PT Informações do produto e instruções / EL Πληροφορίες προϊόντος και οδηγίες / PL Informacje o produkcie oraz instrukcje / CS Informace o výrobku a návod k použití / HU A termék ismertetése és használat / TR Ürün Bilgileri ve Talimatları / NO Produktinformasjon og instruksjoner

Instruments Chirurgicaux Manuels / Manual Surgical Instruments / Strumenti Chirurgici Manuali / Manuelle Chirurgische Instrumente / Instrumentos Quirurgicos Manuales / Handmatige Chirurgische Instrumenten / Manuelle Kirurgiske Instrumenter / Kirurgisten Instrumenttien Kayttoohjeet / Manuella Kirurgiska Instrument / Instrumentos Cirurgicos Manuais / Χειριζόμενα Χειρουργικά Εργαλεία / Ręczne Narzędzia Chirurgiczne / Ručni Chirurgičke Nastroje / Sebeszeti Kezi Műszerek / Manuelle Kirurgiske Instrumenter / Manuel Cerrahi Aletle

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EN English / Product Information and Instructions

DESCRIPTION OF THE DEVICES

Manual Surgical Instruments: They include instrument trays. Physicians should seek relevant pre-clinical education before attempting new surgical procedures, especially those in insufflated cavities.

MATERIALS

The instruments are fabricated of one of the following biocompatible materials: glass, ceramic, titanium nitride, stainless steel, tungsten, thermoset polymers (including silicone), thermoplastic polymers, sterling silver or chromed brass.

INTENDED USE

Manual surgical instruments are intended for use in a wide variety of surgical procedures including general surgery, or specialized open or laparoscopic surgeries

INDICATIONS FOR USE

Manual surgical instruments are designed to perform a specific function, such as cutting, grasping, clamping, dissecting, probing, retracting, draining, aspirating, suturing, or ligating. For use by, or as directed by a surgeon. Instruments should be used only for the purpose for which they are designed. The proper surgical technique for the use of instruments is the responsibility of the surgeon.

CONTRAINDICATIONS

None known.

SAFETY

Before each use, check the outer surface of the parts introduced into the patient to verify that there is no rough surface, no sharp edges or undesirable protrusion which may cause a safety hazard.

WARNINGS

- Do not use these instruments for purposes other than those for which they are intended.
- Do not bend, disassemble or force an instrument; breakage or failure of the instrument could occur, resulting in possible harm to the patient or user.
- Use appropriate personal protection equipment when handling contaminated instruments.
- Do not place the instrument on the patient when it is not being used.
- Do not try to modify the instrument.
- Trocar canulas, if used with the device, should be inspected for sharp edges that may damage the electrical insulation of the instrument. When using a trocar canula, the working end of the device must pass entirely through the trocar canula for the device to function properly and effectively.
- The pathways of the current through conductive elements like metal instruments and endoscopes can cause local burns to the patient, the physician or another member of the care team. Contacting conductive elements with the active cautery area may cause undesired tissue heating and burns.
- To ensure proper functioning of dismantable devices and devices with accessories, check the assembly and the functionality of all elements of the device before use.
- For mouth gags, which include tooth liners, inspect the liners for signs of wear. If the condition of the liners is questionable, replace them before the next procedure.
- Discard instrument after suspected Creutzfeldt-Jakob Disease (CJD) exposure; Integra instruments have not been validated to withstand the chemical and thermal exposure recommended to eradicate prions.

PRECAUTIONS

- Do not use abrasive cleaning products such as hard brushes, etc.
- Extreme care must be taken when handling and cleaning instruments with sharp points or edges.
- We recommend you to use MICROFRANCE® trocars.
- Delicate instruments should be handled with extreme care to prevent damage.

ADVERSE EVENTS

None known.

INSTRUCTION FOR USE

- Reusable instruments must be cleaned and sterilized before the first use and before every new use.
- It is the responsibility of the surgical team to select the appropriate instruments for each use.
- Carefully insert and remove the instruments from the trocars. Close the jaws of the instruments when passing through trocars.
- Check the condition of instruments before and after each case. Remove from use any incomplete or poorly operating instruments.

REPROCESSING INSTRUCTIONS for Reusable Instruments (PER ISO17664)

When applying dry times to Integra instrument trays and their accessories, dry times outside the standard healthcare prevacuum parameters may be required. The current recommended dry times may range from a standard 15 minutes to an extended time of 40 minutes. This is especially important for polymer-based (plastic) instrument trays. The user should have adequate means to verify drying time, such as visual inspection.

| Warnings/ Precautions | To prevent stains, use distilled or demineralized water, and uses a neutral detergent to reprocess these instruments. Do not cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilize, as damage to the instrument finish may occur. | | | | | | | | | | | | | | | | | | | | | | |
|---|--|-------------------|--------------------------------|-------|--------------------|-------------------|----------------|------------|-----------|----------------|-----|-------------|-----------|---------------|--------------------------------|--------|-----------|----------------|----------------------|---------|-----------|---------------|-----|
| INSTRUCTIONS | | | | | | | | | | | | | | | | | | | | | | | |
| 1) Bringing into service of a new instrument | This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. The protection / silicone caps used during shipment must be removed before cleaning and sterilization procedures. | | | | | | | | | | | | | | | | | | | | | | |
| 2) Preparation at the point of use | It is recommended that instruments be reprocessed as soon as is practical following use. Rinse instrument immediately after use in distilled or demineralized water. Keep instrument moist (for example, cover with a wet drape or saturate with tap water – do not use saline) after use and prior to decontamination, cleaning and sterilization to ensure adequate cleaning. Instruments with cleaning ports and/or suction tubes shall be flushed with neutral pH enzymatic detergent, followed by a tap water rinse. | | | | | | | | | | | | | | | | | | | | | | |
| 3) Preparation before cleaning | Fully disassemble modular designed instrument for effective cleaning. Remove any cap covering the cleaning port, if applicable. For instruments with lumens and instruments equipped with cleaning ports, inject cleaning solution through the instrument with an irrigation syringe. When cleaning suction tubes with holes to control suction, place gloved finger over the control hole to flush through the tube. Open any articulated instruments before positioning in the soaking solution. Place instruments so they do not touch each other. Note: Flushing in running water is essential between decontamination and cleaning to prevent any risk of reaction between the two solutions. | | | | | | | | | | | | | | | | | | | | | | |
| 4) Cleaning, disinfection and drying | <p>Cleaning: Automated: All instruments shall be thoroughly rinsed with tap water prior to placement in the automatic washer. Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/disinfector manufacturers. Use alkaline or neutral pH detergent recommended by washer/disinfector or detergent manufacturers. These products have been validated for effective cleaning using an automatic washer/disinfector cycle consisting of the following:</p> <table border="1"> <thead> <tr> <th>Phase</th> <th>Recirculation Time</th> <th>Water temperature</th> <th>Detergent type</th> </tr> </thead> <tbody> <tr> <td>Pre-wash 1</td> <td>3 minutes</td> <td>cold tap water</td> <td>N/A</td> </tr> <tr> <td>Enzyme Wash</td> <td>5 minutes</td> <td>hot tap water</td> <td>neutral pH enzymatic detergent</td> </tr> <tr> <td>Wash 1</td> <td>5 minutes</td> <td>66°C set point</td> <td>neutral pH detergent</td> </tr> <tr> <td>Rinse 1</td> <td>3 minutes</td> <td>hot tap water</td> <td>N/A</td> </tr> </tbody> </table> <p>Thoroughly examine instruments for any residual soil. A rigid MIS (minimally invasive surgery) rack should be used to clean rigid MIS instruments and instruments with lumens to maximize efficiency of cleaning. Suction ports and flush ports should be connected to the MIS Rack to allow a constant flow of solution inside hard-to-reach areas of the instruments.</p> | | | Phase | Recirculation Time | Water temperature | Detergent type | Pre-wash 1 | 3 minutes | cold tap water | N/A | Enzyme Wash | 5 minutes | hot tap water | neutral pH enzymatic detergent | Wash 1 | 5 minutes | 66°C set point | neutral pH detergent | Rinse 1 | 3 minutes | hot tap water | N/A |
| Phase | Recirculation Time | Water temperature | Detergent type | | | | | | | | | | | | | | | | | | | | |
| Pre-wash 1 | 3 minutes | cold tap water | N/A | | | | | | | | | | | | | | | | | | | | |
| Enzyme Wash | 5 minutes | hot tap water | neutral pH enzymatic detergent | | | | | | | | | | | | | | | | | | | | |
| Wash 1 | 5 minutes | 66°C set point | neutral pH detergent | | | | | | | | | | | | | | | | | | | | |
| Rinse 1 | 3 minutes | hot tap water | N/A | | | | | | | | | | | | | | | | | | | | |

| | Cleaning: Manual: Soak in lukewarm (less than 43°C), mild (pH 7.0 – 8.5), enzymatic detergent and deionized water for a minimum of 2 minutes. For instruments with lumens and instruments equipped with cleaning ports, inject cleaning solution through the instrument with an irrigation syringe and allow to soak for a minimum of 2 minutes. Use a soft instrument brush to scrub instruments while submerged in cleaning solution to remove organic matter. Rinse with deionized water, then clean ultrasonically in a lukewarm (less than 43°C), mild (pH 7.0 – 8.5) detergent and deionized water for 10 minutes. Rinse thoroughly with deionized water, utilizing a syringe to thoroughly rinse cleaning solution from lumens and cleaning ports. Clean in this manner until no visible soil remains on the instrument. Dry with compressed air or wipe dry with a lint-free cloth. Examine instruments for any staining or deterioration, remove from use as appropriate. Note: When using an ultrasonic cleaner or a spray washing machine, follow the manufacturer's recommendations, particularly with regard to articulated instruments and positioning of instruments. Note: Additional cleaning methods may be warranted, including presoaking in 3% hydrogen peroxide. | | | | | | | | | | | | | | | | | | | | | | | | |
|---|---|---------|---|------------------|---|--|--|-------|---------|---------|---------|------------------|----------------|-------------|--------|--------|--------|--------|--------|------|--------|--------|-------|--------|-------|
| 5) Maintenance, inspection and testing | Check the cleanliness and operation of the instrument. Clean again if debris is present and remove from use any damaged instrument. Inspect components for any damage. If damage is observed, do not use the instrument until it is repaired. Following cleaning, lightly lubricate instruments with movable parts. Use a lubricant intended for sterilizable instruments such as a water-soluble instrument milk. Do not use silicone spray. | | | | | | | | | | | | | | | | | | | | | | | | |
| 6) Packaging | A standard sterilization wrap may be used. In the US, an FDA cleared surgical wrap is required. In sets: Instruments may be loaded into dedicated instrument trays or general purpose sterilization trays. Ensure that cutting edges are protected. Wrap trays using appropriate method. | | | | | | | | | | | | | | | | | | | | | | | | |
| 7) Sterilization (Temperatures are minimum required; times are minimum required) | The sterilization parameters given below should be used for devices that are fully disassembled when disassembly is possible. During sterilization, do not mix instruments of different metals, to avoid electrolyte oxidation. Close instruments with catches and racks on the first notch. Arrange the instruments in sterilization containers with perforations on the top and bottom and on supports such as those used in microsurgery. Follow the appropriate cycle listed in the table below. All steam cycles have been validated in the wrapped configuration and can be sterilized wrapped or unwrapped. These devices have only been validated for steam sterilization methods. | | | | | | | | | | | | | | | | | | | | | | | | |
| | <table border="1"> <thead> <tr> <th></th> <th colspan="2">Instruments Only</th> <th colspan="3">Instruments Only or Instruments in Instrument Trays</th> </tr> <tr> <th>Cycle</th> <th>Gravity</th> <th>Gravity</th> <th>Pre-vac</th> <th>Pre-vac (FR/WHO)</th> <th>Pre-vac (U.K.)</th> </tr> </thead> <tbody> <tr> <td>Temperature</td> <td>121 °C</td> <td>132 °C</td> <td>132 °C</td> <td>134 °C</td> <td>134 °C</td> </tr> <tr> <td>Time</td> <td>30 min</td> <td>10 min</td> <td>4 min</td> <td>18 min</td> <td>3 min</td> </tr> </tbody> </table> | | Instruments Only | | Instruments Only or Instruments in Instrument Trays | | | Cycle | Gravity | Gravity | Pre-vac | Pre-vac (FR/WHO) | Pre-vac (U.K.) | Temperature | 121 °C | 132 °C | 132 °C | 134 °C | 134 °C | Time | 30 min | 10 min | 4 min | 18 min | 3 min |
| | Instruments Only | | Instruments Only or Instruments in Instrument Trays | | | | | | | | | | | | | | | | | | | | | | |
| Cycle | Gravity | Gravity | Pre-vac | Pre-vac (FR/WHO) | Pre-vac (U.K.) | | | | | | | | | | | | | | | | | | | | |
| Temperature | 121 °C | 132 °C | 132 °C | 134 °C | 134 °C | | | | | | | | | | | | | | | | | | | | |
| Time | 30 min | 10 min | 4 min | 18 min | 3 min | | | | | | | | | | | | | | | | | | | | |
| | Drying: Instruments Only and Aluminum/Stainless Steel Instrument Trays : 15 to 30 minutes and visual inspection for dryness. Polymer-based (plastic) Instrument Trays: Minimum drying time 40 minutes and visual inspection for dryness | | | | | | | | | | | | | | | | | | | | | | | | |
| 8) Storage | Store instruments in a clean, dry area. | | | | | | | | | | | | | | | | | | | | | | | | |
| 9) Disposal | Instruments must be disposed of in accordance with applicable local regulations. | | | | | | | | | | | | | | | | | | | | | | | | |

The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. They are not applicable to single use devices or single use accessories, which must be destroyed after use in accordance with applicable local regulations. It remains the responsibility of the processor to ensure that reprocessing is performed using validated equipment to achieve the desired result. This normally requires validation and routine monitoring of the process. Some devices have specific assembly instructions. In this case, refer to the assembly insert provided with the device for additional instructions. All validations performed per AAMI TIR12, (Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers).

RETURNS

Contact Integra Customer Service to obtain a Return Goods Authorization number (RGA #) prior to shipping the product back to Integra. The RGA # should be prominently displayed on the outside of the return packaging and included on all paperwork enclosed with the return. All product returned should be decontaminated and cleaned, and all products should be safely packed in protective wrapping. Please have the original invoice number or purchase order number available to assist in verifying warranty information.

LIMITED WARRANTY

A. This LIMITED WARRANTY provides assurance for the customer who purchases a Integra Product (hereinafter the "Product") that should the Product fail to function to Integra's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product), Integra will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Integra or from its affiliate or its authorized distributor or representative.

B. To qualify for this LIMITED WARRANTY, the following conditions must be met:

- (1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
- (2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
- (3) Integra must be notified in writing within thirty (30) days following discovery of a defect.
- (4) The Product must be returned to Integra within thirty (30) days of Integra receiving notice as provided for in (3) above.
- (5) Upon examination of the Product by Integra, Integra shall have determined that: (i) the Product was not repaired or altered by anyone other than Integra or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Integra be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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