

Integra®

CUSA® NXT
Ultrasonic Tissue Ablation System

Limit uncertainty by
simplifying operation and
maximizing performance
with NXT Digital Architecture™



INTEGRA®
LIMIT UNCERTAINTY

CUSA® NXT Ultrasonic Tissue Ablation System

The CUSA NXT is the newest member of the CUSA® family. The CUSA NXT system features revolutionary NXT Digital Architecture™ and a touch-screen user interface.

The CUSA NXT supports the current selection of the CUSA Selector handpieces. When used with the CUSA Selector MRI handpieces, the system is fully compatible for use in suites with intraoperative MRI equipment.

The CUSA NXT is a compact and portable unit. The Console can be connected to the wall suction in the operating room, or used with the CUSA NXT Service Module for optimized aspiration performance.




CUSA® NXT System

Catalog#	Description	Quantity
CUSANXT	CUSA NXT Console, Footswitch, Power Cord, Operator's Manual	Each 1
CUSANXT1	CUSA NXT Service Module	Each 1

CUSA® Selector® Neuro Handpieces

24 kHz Neuro Long Handpiece

Catalog#	Sterile	Description	Quantity
1523000M5	No	24 kHz Neuro Long Handpiece	Each 1
1523225	Yes	24 kHz Neuro Long Sterile Tip Set and Flue	Box 5
1523067	No	24 kHz Neuro Non-Sterile Tip Set and Flue	Box 5
1523049	No	24 kHz Neuro Long Angled Shroud	Each 1
1523045	No	24 kHz Neuro Long Clear Extension Flue	Box 5
1523052	No	24 kHz Neuro Long Angled Extension	Each 1




1-Transducer assembly, 2-Angled extension,
3-24 kHz Neuro Long Tip Set (single use) sterile, 4-Angled shroud, 5-Extension flue

24 kHz Neuro Short Handpiece

Catalog#	Sterile	Description	Quantity
1523000M7	No	24 kHz Neuro Short Handpiece	Each 1
1528020M7	No	24 kHz Neuro Short Handpiece for MRI	Each 1
1523215	Yes	24 kHz Neuro Short Sterile Tip Set and Flue	Box 5
1523067	No	24 kHz Neuro Short Non-Sterile Tip Set and Flue	Box 5
1523216	Yes	24 kHz Selector SaberTip Set and Flue	Box 5
NXT3218*	Yes	24 kHz CUSA NXT Superior Forward BoneTip	Box 5
NXT3219*	Yes	24 kHz CUSA NXT Superior Reverse BoneTip	Box 5
1523049	No	24 kHz Neuro Short Angled Shroud	Each 1
1523053	No	24 kHz Neuro Short Angled Extension	Each 1

* Compatible ONLY with the CUSA NXT System




1-Transducer assembly, 2-Angled extension,
3-24 kHz Neuro Short Tip Set (single use) sterile, 4-Angled shroud

CUSA® Selector® Neuro Handpieces

35 kHz Neuro Handpiece


Catalog#	Sterile	Description	Quantity
1523000M9	No	35 kHz Neuro Handpiece	Each 1
1528020M9	No	35 kHz Neuro Handpiece for MRI	Each 1
1523235	Yes	35 kHz Neuro Sterile Tip Set and Flue	Box 5
1523078	No	35 kHz Neuro Non-Sterile Tip Set and Flue	Box 5
1523243	Yes	35 kHz CUSA NXT Extended Length Tip	Box 5
1523060	No	35 kHz Neuro Angled Shroud	Each 1
1523058	No	35 kHz Neuro Angled Extension	Each 1



1-Transducer assembly, 2-35 kHz Neuro Tip Set (single use) sterile,
3-Angled extension, 4-Angled shroud

24 kHz Neuro Micro Handpiece

Catalog#	Sterile	Description	Quantity
1529000M2	No	24 kHz Micro Handpiece	Each 1
1532020M2	No	24 kHz Micro Handpiece for MRI	Each 1
1529165	Yes	24 kHz Micro Sterile Tip Set and Flue	Box 5
1529066	No	24 kHz Micro Non-Sterile Tip Set and Flue	Box 5
1529024	No	24 kHz Micro Angled Shroud	Each 1
1529043	No	24 kHz Micro Angled Extension	Each 1




1-Transducer assembly, 2-24 kHz Micro Tip Set (single use) sterile,
3-Angled extension shroud, 4-Angled shroud

CUSA® Selector® Standard Handpieces

24 kHz Straight Handpiece


Catalog#	Sterile	Description	Quantity
1523000M1	No	24 kHz Straight Handpiece	Each 1
1523013	No	24 kHz Non-Sterile Standard Tip	Each 1
1523042	No	24 kHz Flue Pack	Box 20
1523017	No	24 kHz Straight Shroud	Each 1



1-Transducer assembly, 2-Flue, disposable,
3-Standard shroud, 4-Standard tip

24 kHz Angled Handpiece

Catalog#	Sterile	Description	Quantity
1523000M2	No	24 kHz Angled Handpiece	Each 1
1523015	No	24 kHz Non-Sterile Standard Tip	Each 1
1523042	No	24 kHz Flue Pack	Box 20
1523049	No	24 kHz Angled Shroud	Each 1
1523033	No	24 kHz Angled Extension	Each 1
1523035	No	24 kHz Extension Shroud	Each 1




1-Transducer assembly, 2-Angled extension, 3-Standard tip,
4-Angled shroud, 5-Extension shroud, 6-Flue, disposable

CUSA® Selector® Standard Handpieces

35 kHz Straight Handpiece


Catalog#	Sterile	Description	Quantity
1523000M3	No	35 kHz Straight Handpiece	Each 1
1523014	No	35 kHz Non-Sterile Standard Tip	Each 1
1523043	No	35 kHz Straight Clear Flue Pack	Box 20
1523018	No	35 kHz Straight Shroud	Each 1



1-Transducer assembly, 2-Flue, disposable, 3-Standard shroud, 4-Standard tip

24 kHz Laparoscopic Handpiece

Catalog#	Sterile	Description	Quantity
1523000M6	No	24 kHz Laparoscopic Handpiece	Each 1
1523029	No	24 kHz Non-Sterile Laparoscopic Tip	Each 1
1523032	No	24 kHz Laparoscopic Flue Pack	Box 20
1523028	No	24 kHz Step Amplifier	Each 1
1523030	No	24 kHz Laparoscopic Shroud	Each 1
1523031	No	24 kHz Laparoscopic Shroud Extension	Each 1



1-Transducer assembly, 2-Flue, reusable, 3-Shroud extension, 4-Tip, 5-Shroud, 6-Step amplifier, 7-Tubing clamp

CUSA® Selector® and CUSA® NXT Accessories

Catalog#	Sterile	Description	Quantity
1517079	Yes	CUSA Selector Sterile Tubing Kit, Single Use	Box 10
1529036	Yes	CUSA Selector Sterile Micro Tubing Kit Extension For Use with Microsurgical Handpiece, Single Use	Box 10
1532024	Yes	CUSA Selector Sterile MRI Tubing Kit Extension For Use with MRI Handpieces, Single Use	Box 10
NXT00005	No	CUSA NXT Contamination Guard	Box 2
07289	No	CUSA Selector Handpiece Sterilization Tray	Each 1
1523101	No	CUSA Selector Handpiece Aspiration Port Cleaning Brush	Each 1

CUSA® Selector® and CUSA® NXT Replacement Parts

Catalog#	Sterile	Description	Quantity
1523072M3	No	Wrench Set Neuro Handpiece	Each 1
1529056	No	Wrench Set Micro Handpiece	Each 1
1523072M1	No	Wrench Set Standard Handpiece	Each 1
1523072M2	No	Wrench Set Laparoscopic Handpiece	Each 1
NXT10001	No	CUSA NXT Footswitch	Each 1
1520000M1	No	CUSA Selector Footswitch	Each 1
1528030	No	CUSA Selector MRI Footswitch	Each 1
NXT90001		CUSA NXT Operator's Manual	Each 1
OM1530000M3		CUSA Selector Operator's Manual	Each 1
OM1532000M1		CUSA Selector MRI Operator's Manual	Each 1



Integra®

CUSA® NXT Ultrasonic Tissue Ablation System

For more information or to place an order, please contact:
Integra ▪ 311 Enterprise Drive, Plainsboro, NJ 08536
800-997-4868 USA ▪ 609-936-5400 outside USA
888-980-7742 fax
integralife.com

