

Integra™

MLX Surgical Headlight System

Limit uncertainty with
maximum light to visualize
critical anatomy



INTEGRA™
LIMIT UNCERTAINTY

The whitest and brightest Xenon light from Integra

15% brighter

True 1000 hour warranty on Xenon lamp (not pro-rated)

New Smart Fan technology runs quieter and cooler

- Temperature-sensing feature adjusts Smart Fan speed to compensate for operating room air temperature
- Cooler running extends lamp life

MLX lamp stays brighter longer than standard 300W lamps*

1 – Stand-by Button


- Eliminates “stray” light when headlight is unplugged
- Saves lamp hours by automatically shutting off lamp after 15 minutes in “stand by” mode

2 – LCD Screen

- Displays lamp and system hours as well as system diagnostics

3 – Membrane Activation Switch

- Digitally controls light intensity in accurate 5% increments
- Activates display of system status



* Data on file at Integra

UltraLite® Pro – Maximum Comfort

1 – Occipital Basket


- Balances and distributes the weight of the headlight and reduces neck strain

2 – Cranial Support Pad

- Stabilizes the headlight

3 – Adjustable Headband

- Provides perfect sizing without temporal compression



Maximum Headlight Performance

Superior optics

- Optimized light output*
- Superior color fidelity

Enhanced Design

- New, lighter design
- 3-Point Support System for comfortable stability
- Extended linkage for better vertical alignment with surgical loupes



UltraLite® Premium Headlight Cable

- Lightweight and flexible 3mm cable
- Bifurcated design allows even weight distribution and balance
- Translucent outer sheath provides visualization of fiber integrity
- Lifetime Warranty Against Browning

UltraLite® Pro Headlight and Accessories

Catalog Number	Description
AX2100BIF	UltraLite Pro headlight with 9ft (275cm) premium bifurcated cable, headlight module and gown clips
001370LX	UltraLite Pro replacement headband
001380LX	Headlight module, black
001388LX9	UltraLite Pro premium bifurcated headlight cable, 9ft (275cm)
001387LX	UltraLite Pro headlight replacement pad kit

MLX Light Source and Accessories

Catalog Number	Description
00MLX	MLX 300W Xenon Light Source with 20' US power cord
00MLXEU	MLX 300W Xenon Light Source with 6.1m EU power cord
00MLXAU	MLX 300W Xenon Light Source with 6.1m Australian power cord
00MLXUK	MLX 300W Xenon Light Source with 6.1m UK power cord
001320LX	Lamp Module Replacement MLX
001337	Floorstand for MLX 300W Xenon Light Source

For more information or to place an order, please contact:
Integra ■ 311 Enterprise Drive, Plainsboro, NJ 08536
USA and Canada: 800.431.1123 ■ 914.592.9050 (Outside USA)
914.592.8056 (Fax)
integralife.com